



415 S Airpark Road
Cottonwood, AZ 86326, USA

Phone: (800) 733-0266
Fax: (928) 649-2306



F25/26 Moulton Park Business Centre,
Redhouse Road
Northampton, NN3 6AQ, UK

Phone: 01604 497516
Fax: 01604 497501

NUT/RB

MICROSLIDE® TECHNICAL DOCUMENT



Index

USE	Page 3
APPLICATION	Page 3
PADDLE AGARS	Page 3
STORAGE/EXPIRATION	Page 3
AGAR VERIFICATION	Page 4
SAMPLING	Page 4
INCUBATION	Page 5
COLONY MEASURING	Page 5
ENUMERATION	Page 5
DISPOSAL	Page 6
IDENTIFICATION	Page 6
GLOSSARY	Page 11

NUT/RB

CODE: M-NUT/RB

USE

Isolation and differentiation of Gram (-) enteric bacilli. **(NUT)** Selective enumeration and cultivation of yeasts, molds, and Actinomycetes from food and other surfaces **(RB)**.

APPLICATION

In total coliform testing (TCC), the coliform organisms tested for include: total coliform, fecal coliform, and *E. coli* (*Escherichia coli*). Detection of fecal coliforms (a subset of total coliforms) or *Escherichia coli* (a subset of fecal coliforms) can indicate the potential presence of waterborne pathogens associated with fecal contamination¹. Rose Bengal Agar is recommended in *Standard Methods* for the enumeration of yeasts and molds from food and water.

PADDLE AGARS



Side 1: Nutrient-TTC Agar (NUT) – (Color: Yellow) General purpose (relatively non-selective) medium, which will support the growth of a wide variety of organisms. Suitable for cultivation of both aerobes and anaerobes. Aerobic coliform bacteria can be detected by their ability to reduce the TTC dye to a red-colored formozan dye. Bacterial colonies appear as red dots on an otherwise yellow medium.

Note: The Nutrient-TTC agar color is normally light yellow when the agar is cast. After testing, during the incubation phase, the agar may change to a light green color. This color change is a result of either a microbial-induced or chemically-induced pH change in the media. This color change alone does not indicate the presence of microorganisms. Development of red spots or other growth on the agar are an indication of microorganisms.



Side 2: Rose Bengal Agar (RB) – (Color: Pink) Selective medium for the enumeration of yeasts and molds.

***Note:** Side 1 of each paddle is marked with an indented laser line.

STORAGE / EXPIRATION

Microslides[®] should be stored tightly sealed (unopened) in a cool, dry location at room temperature (18 - 25°C; 65 - 77°F). Temperature fluctuations may result in condensation settling at the bottom of the vial, although this does not affect culture properties, it could reduce the shelf-life or cause the agar to separate from the plastic paddle support. Refer to 'Best Before End date' (SEE: BBE stamped on vial).

Avoid sudden temperature changes. Shield from direct sunlight. Do not allow paddles to freeze. Do not store in a refrigerator (~44°F / 10°C) or at temperatures exceeding 80°F; 27°C. Refrigeration may result in

¹ United States Pharmacopeial Convention. 2007. The United States pharmacopeia, 31st ed., Amended Chapters 61, 62, 111. The United States Pharmacopeial Convention, Rockville, MD.

water condensation. Discard if paddle agar appears oxidized (darkened from expected color) or if contaminants appear. Expiry applies to medium in its intact container when stored as directed.

AGAR VERIFICATION

These agars have been verified by [EMSL Analytical, Inc.](#) using *E. coli* and *E. faecalis* (NUT) and *P. commune* and *C. albicans* (RB) cultures. Documentation available upon request.

SAMPLING

SURFACE Sampling Protocol

1. Remove the paddle from the vial. Do not touch the agar surfaces.
2. To assure an accurate area recovery, contact the paddle to 20²cm of the surface by contacting the surface twice in separate 10²cm areas.
3. Replace paddle in vial.
4. Incubate.

LIQUID Sampling Protocol

DIRECT IMMERSION PROTOCOL – low viscous liquids

1. Mix liquid test sample.
2. Remove the paddle from the vial. Do not touch the agar surfaces.
3. When taking the sample:
 - a. Pour 40mL of the sample into the vial (to the printed horizontal fill line; see right). Dip the paddle into the 40mL volume liquid in the vial. Maintain a contact time of at least 15 seconds (30 seconds optimal). Both agar surfaces must be completely contacted.
 - b. Or dip the paddle into the sample directly. Maintain a contact time of at least 15 seconds (30 seconds optimal). Both agar surfaces must be completely contacted.
4. Allow excess fluid to drain off both paddle agar surfaces.
5. Replace paddle in vial.
6. Incubate.



SPREAD Protocol – high viscous liquids

1. Mix liquid test sample.
2. Remove paddle from vial. Do not touch the agar surfaces.
3. Holding the contact agar surface on a horizontal plane, deposit volume as a single drop approximately 1cm from the handle boundary (Figure 1).
4. Position a sterile glass rod on the "handle" side of the drop and bring it into contact with the drop creating a meniscus. Drag the glass tube over the paddle agar surface.
5. Replace paddle in vial.
6. Incubate.




Figure 1


INCUBATION

Incubation of Paddle Growth	Incubation Temperature	Examine at:
Yeast / Mold	25 to 30°C	48 hours up to 120 hours (5 days)
Yeast / Mold	Room Temperature	Up to 7 days
Total Coliform / Bacteria	35 ± 2°C	24 to 48 hours
Total Coliform / Bacteria	Room Temperature	Up to 5 days

Note: Incubation of bacteria after 48 hours may produce confluent growth making enumeration more difficult.

COLONY MEASURING

Each Microslide® paddle has molded media attachment points that are 4mm in length (point-to-point). This feature provides a useful guidepost to estimating nearby colony size.




ENUMERATION


Bacteria CFU/mL




10³ cfu/mL
(1,000)
(Light)




10⁴ cfu/mL
(10,000)



10⁵ cfu/mL
(100,000)
(Moderate)



10⁶ cfu/mL
(1,000,000)




10⁷ cfu/mL
(10,000,000)
(Heavy)

Note: Estimation of lower counts is possible, but statistically difficult to justify. Use Light, Moderate and Heavy for Mold growth and surface testing.







DISPOSAL

Make a 1:9 dilution of household bleach (5.25% sodium hypochlorite solution). Twist and remove Microslide® paddle from vial. Fill vial with 40mL diluted hypochlorite solution (to fill-line). Allow 15-minute contact time. Discard bleach solution. Replace paddle in vial and dispose. Alternatively, loosen cap and microwave for 30 seconds, autoclave, or incinerate.

IDENTIFICATION

Organism	Nutrient-TTC (NUT)	Rose Bengal (RB)
<i>Actinomyces bovis</i>	Growth: + Colony: Opaque/tan-grey, CVEG, 1-3mm	Growth: ++ Colony: Opaque/tan-grey, CVEG, 1-3mm
<i>Alternaria spp.</i>	Growth: + Colony: Downy to woolly; flat, grayish, short, aerial hyphae, later becomes greenish black or olive-brown with a light border, 3-9cm	Growth: ++ Colony: Suede-like to woolly, initially white to yellow-orange, becoming black to olive-green or grayish, or grayish-green, umbonate with lighter center area, condensation (rings), fast-growing, 3-9cm+ (confluent growth)
<i>Aspergillus niger</i>	 Growth: +++ Colony: Granular, jet black conidia with yellow/gray hyphae, 3-5+cm	 Growth: +++ Colony: Woolly and/or felt-like, forms a carpet, initially white later with jet black fruiting bodies (sporangia), fast-growing (4.5cm in 4 days), 3-9cm+ (confluent growth)
<i>Aspergillus flavus</i>	Growth: + Colony: Granular to woolly, yellow, yellow-green, or yellow-brown, 3-9cm	Growth: +++ Colony: Granular to woolly, yellow, yellow-green, or yellow-brown, 3-9cm+ (confluent growth)
<i>Aspergillus fumigatus</i>	Growth: + Colony: Granular to cottony, blue-green, green-grey, or green-brown, 3-9cm	Growth: +++ Colony: Felt-like, forms a carpet, initially white to green or blue-green fruiting bodies, 3-9cm+ (confluent growth)
<i>Aspergillus terreus</i>	Growth: + Colony: Granular, radially rugose (wrinkled), cinnamon buff/brown, 3-9cm	Growth: +++

<p><i>Bacillus spp.</i></p>	 <p>Growth: +++ Colony: Green with dark green center</p>	<p>Colony: Granular, radially rugose (wrinkled), cinnamon buff/brown, 3-9cm+ (confluent growth)</p>  <p>Growth: ++ Colony: Translucent to pink, circular to irregular, flat to raised, entire, 2-5mm</p>
<p><i>Botrytis spp.</i></p>	<p>Growth: + Colony: Woolly, white/grey/brown pigment, 3-9cm</p>	<p>Growth: +++ Colony: Woolly, white/grey/brown pigment, 3-9cm</p>
<p><i>Candida albicans</i></p>	 <p>Growth: +++ Colony: Cream, CVEG, 1-2mm</p>	 <p>Growth: +++ Colony: White to pink, circular, convex, dull, entire, 0.1-0.5mm</p>
<p><i>Chaetomium spp.</i></p>	<p>PARTIAL TO COMPLETE INHIBITION</p>	<p>Growth: +++ Colony: Suede-like to Woolly, initially white, later globular (roundish) gray or olive areas / structures (perithecia) looking like cockleburs, 3-5cm+ (confluent growth)</p>
<p><i>Cladosporium spp.</i></p>	 <p>Growth: + Colony: Granular to woolly (velvety), olive-brown to black/brown, sometimes grey on a dark base, 2-5+++cm</p>	 <p>Growth: + Colony: Suede-like to woolly, often becoming powdery due to the production of abundant conidia, forms a carpet, white turning olive-brown, buff, or brown, slow-growing, 3-9cm+ (confluent growth)</p>
<p><i>Epicoccum spp.</i></p>	<p>Growth: + Colony: Woolly, cottony, felty, yellow/orange/red, 3-5cm</p>	<p>Growth: +++ Colony: Woolly, cottony, felty, yellow/orange/red, 3-5cm</p>

<i>E. coli</i>	 <p>Growth: +++ Colony: Yellow/Orange/Red, CVEG, 0.5-1.0mm</p>	INHIBITED
<i>Enterobacter aerogenes</i>	 <p>Growth: +++ Colony: Red, CVEG, 2-4mm</p>	Growth: ++ Colony: Pink to red, CVEG, 2-4mm
<i>Enterococcus spp.</i> <i>Fusarium spp.</i>	 <p>Growth: + Colony: Wooly, flat, sometimes mucous-like</p>	INHIBITED  Growth: +++ Colony: Wooly, initially white, later with yellow, pink, red, purple or pale brown coloring, fast-growing, 3-9cm+ (confluent growth)
<i>Klebsiella spp.</i>	 <p>Growth: +++ Colony: Amber/Red, spreading, 4-5mm</p>	INHIBITED
<i>Microsporium spp.</i>	 <p>Growth: + Colony: Glabrous (smooth), downy, wooly, powdery, white at first, later becoming grayish-yellow to blue-green with age, wrinkled with age, 1-9+cm</p>	Growth: + Colony: Glabrous (smooth), downy, wooly, powdery, white at first, later becoming grayish-yellow to blue-green with age, wrinkled with age, 1-9+cm

<i>Mucor spp.</i>	 <p>Growth: + Colony: Woolly, fluffy (like cotton candy), white at first, later becoming gray/yellow to blue-green with age, 2-5++cm</p>	 <p>Growth: + Colony: Woolly, initially white, then white-yellow to various shades of gray to green with lollipop fruiting bodies (sporangia), fast-growing, 3-9cm+ (confluent growth)</p>
<i>Penicillium chrysogenum (notatum)</i>	 <p>Growth: ++ Colony: Granular, velvety/powdery, flat, initially white, then various shades of green-blue, green, or yellow-green, 3-5cm</p>	 <p>Growth: ++ Colony: Granular, velvety/powdery, flat, initially white, then various shades of green-blue, green, or yellow-green, 3-9cm+ (confluent growth)</p>
<i>Penicillium roqueforti</i>	 <p>Growth: + Colony: Granular, dull, green in colour, arachnoid (with many spider web-like fibers) colony margins, 0.5-1.0cm</p>	 <p>Growth: ++ Colony: Granular, velvet-like, flat, initially white then various shades of green, blue-green pigment, 3-9cm+ (confluent growth)</p>
<i>Penicillium digittum</i>	<p>Growth: + Colony: Woolly, fluffy (like cotton candy), white at first, later becoming green with age, 3-9cm</p>	<p>Growth: +++ Colony: Suede-like, woolly, initially white, then various shades of olive green, 3-9cm+ (confluent growth)</p>
<i>Pithomyces spp.</i>	<p>Growth: + Colony: Powdery, pale/dark grey or brown pigment, 2-9++cm</p>	<p>Growth: +++ Colony: Powdery, pale/dark grey or brown to olive green pigment, lighter outer ring with center bullseye, 2-9cm+ (confluent growth)</p>
<i>Proteus spp.</i>	INHIBITED	INHIBITED

Pseudomonas aeruginosa



INHIBITED

Growth: +++
Colony: Red, irregular, spreading to confluent, 2-4mm

Pseudomonas fluorescens



Growth: +++
Colony: Clear/colorless with grey/dark center, translucent edges, irregular/spreading to confluent, 2-4mm



Growth: +
Colony: Translucent, pinkish, or amber, irregular, raised, undulate, 2-4mm+

Rhizopus spp.



Growth: +++
Colony: Cottony, white to black/grey (black fruiting bodies), 2-9++cm



Growth: +++
Colony: Dense, cottony growth, initially white, turning to gray with black fruiting bodies (sporangia), fast-growing, 3-9cm+ (confluent growth)

Saccharomyce cerevisiae



Growth: ++
Colony: Creamy white to tan, spreading, circular, entire, raised to convex, glistening surface, 5-8mm









Growth: +++
Colony: Translucent to white or cream, CVEG (may be dull), 0.1-0.5mm (punctiform)

Salmonella typhimurium

Growth: +++
Colony: Purple/pink, FED, 0.5-1.0mm

INHIBITED

<i>Salmonella epidermidis</i>		INHIBITED
<i>Serratia spp.</i> <i>Shigella spp.</i>		PARTIAL TO COMPLETE INHIBITION INHIBITED
<i>Staphylococcus aureus</i>		INHIBITED
<i>Streptococcus spp.</i>		INHIBITED
<i>Torula spp.</i> <i>Trichoderma spp.</i>		 Growth: +++ Colony: White, opaque, viscous, CVEG, 0.1-0.5mm (punctiform) Growth: ++

For *in vitro* diagnostic use only. This product should be used only by adequately trained personnel with knowledge of microbiological techniques in the laboratory. © Precision Laboratories, Inc. All rights reserved.

<i>Trichophyton</i> <i>spp.</i>	Colony: Cottony, white, later scattered green or yellow-green patches (rings), 2-9++cm	Colony: Cottony, white, later scattered green or yellow-green patches (rings), 3-9cm+ (confluent growth)
	Growth: + Colony: Wooly with indented borders, white to brown/tan pigment, 2-9++cm	Growth: ++ Colony: Wooly, initially white with brownish/tan pigmentation, outer darker ring, 3-9cm+
<i>Gram (+)</i> <i>Bacteria</i>	PARTIAL TO COMPLETE INHIBITION	

GLOSSARY

CVEG.....Convex, Entire, Glossy

FED.....Full, Entire, Dull

Gram.....Gram reaction