





Precision

LABORATORIES





NUTRIENT TTC-ROSE BENGAL DIPSLIDE

SKU: M-NUT-RB

The Nutrient TTC-Rose Bengal Dipslide is often used in testing personal care products/cosmetics, as well as food and water.

Nutrient-TTC Agar is a general purpose (relatively non-selective) medium, which will support the growth of a wide...

[Read More](#)



TRYPTIC SOY-ROSE BENGAL DIPSLIDE

SKU: M-TSA-RB

The Tryptic Soy-Rose Bengal Dipslide is a general purpose paddle often used in educational settings. Tryptic Soy Agar is an enriched media, suitable to support fastidious heterotrophs and to facilitate vigorous growth of aerobic and anaerobic...

[Read More](#)



PH 2-9 TEST STRIPS, 3 PAD

SKU: PH2090-3

The pH 2-9 test strips contain three pads for testing accuracy. This strip measures in steps of 0.5 pH units, and gives results in just seconds. pH 2-9 test strips measures a wide range of both acidic and alkaline solutions. This strip can be...

[Read More](#)



PH 7-14 TEST STRIP, 3 PAD

SKU: PH0714-3

The pH 7-14 test strip is a base pH strip containing three pads for testing accuracy. The pH 7-14 test strip measures in increments of 0.5, and is used to test neutral to very alkaline solutions. Please see notes below regarding how to read the...

[Read More](#)

