





# Precision

LABORATORIES





## SULFITE TEST STRIP

**SKU:** SUL

The Sulfite Test Strip is used to detect Sulfite ions in water-based solutions, measuring in increments of 10, 50, 100, 250, and 500ppm.

The Sulfite test strip will detect free sulfite. Detection of combined or total (free + combined) sulfite...

[Read More](#)



## WINE PH TEST STRIPS 2.8-4.4

**SKU:** PH2844

Wine pH test strips are used to test the pH of the juice prior to fermentation, but can also be used to test the pH of kombucha tea.

[Read More](#)



## BEER PH TEST STRIP 4.6-6.2

**SKU:** PH4662

The narrow-range Beer pH test strip, pH 4.6-6.2, is used for testing beer mash. The pH of Yeast is best at 5.2-5.5.

[Read More](#)

