




# Precision

LABORATORIES





## WOOD SPLINTS

**SKU:** WDS

Wood splints are used for mixing and applying in a laboratory setting.

[Read More](#)

## LENS PAPER, BIBULOUS PAPER

**SKU:** LENS, BIB, BLO

The Lens Book consists of 50 sheets of premium grade white tissue paper. The paper is suitable for general cleaning of glass and plastic lens. It is recommended for cleaning microscope eyepieces and objectives in a classroom setting. This is NOT...

[Read More](#)

