


Precision

LABORATORIES





NUTRIENT TTC DIPSLIDE

SKU: M-NUT-NUT

Nutrient TTC Agar is a general purpose (relatively non-selective) medium, which will support the growth of a wide variety of organisms. Suitable for cultivation of both aerobes and anaerobes. Aerobic coliform bacteria can be detected by their...

[Read More](#)



NUTRIENT TTC-ROSE BENGAL DIPSLIDE

SKU: M-NUT-RB

The Nutrient TTC-Rose Bengal Dipslide is often used in testing personal care products/cosmetics, as well as food and water.

Nutrient-TTC Agar is a general purpose (relatively non-selective) medium, which will support the growth of a wide...

[Read More](#)



NUTRIENT TTC- MACCONKEY DIPSLIDE

SKU: M-NUT-MAC

Nutrient-TTC Agar is a general purpose (relatively non-selective) medium, which will support the growth of a wide variety of organisms. Suitable for cultivation of both aerobes and anaerobes. Aerobic coliform bacteria can be detected by their...

[Read More](#)



NUTRIENT TTC-R2A DIPSLIDE

SKU: M-NUT-R2A

Nutrient-TTC Agar is a general purpose (relatively non-selective) medium, which will support the growth of a wide variety of organisms. Suitable for cultivation of both aerobes and anaerobes. Aerobic coliform bacteria can be detected by their...

[Read More](#)



NUTRIENT TTC-MALT DIPSLIDE

SKU: M-NUT-MALT

The weekly monitoring of bacteria levels in cooling water has been recommended by many legislators and professional authorities as a visual performance indicator to both system and treatment regime. This allows the user to gauge how effective a...

[Read More](#)



TRYPTIC SOY-ROSE BENGAL DIPSLIDE

SKU: M-TSA-RB

The Tryptic Soy-Rose Bengal Dipslide is a general purpose paddle often used in educational settings. Tryptic Soy Agar is an enriched media, suitable to support fastidious heterotrophs and to facilitate vigorous growth of aerobic and anaerobic...

[Read More](#)



RANGER DIPSLIDE INCUBATOR

SKU: INC-R45

Portable incubation of up to 45 dipslides. The Ranger Dipslide Incubator is designed for life in the field or in vehicles with its robust case resistant to oil and chemical environments. It can even be fit to an auxiliary battery so the incubator...

[Read More](#)



CHALLENGER DIPSLIDE INCUBATOR

SKU: INC-C25

Portable incubation of up to 25 dipslides. The Challenger Dipslide Incubator case is easy to wipe clean and robust enough for most hostile work environments. The ergonomic location of the controls under the lid eliminates the risk of accidental...

[Read More](#)



DISCOVERY DIPSLIDE INCUBATOR

SKU: INC-D8

Portable incubation of up to 8 dipslides. The Discovery Dipslide incubator is designed for life in the field or vehicles. The Temperature Control Technology is a fully solid state with no moving parts, and element switching is transistorized and...

[Read More](#)



GNAT DIPSLIDE INCUBATOR

SKU: INC-G2

Portable incubation of 2 dipslides. The Gnat Dipslide Incubator is the easiest to operate. Just plug it in and within 30 minutes, it will heat up to the pre-set temperature and incubate until you turn it off. There are no additional buttons or...

[Read More](#)

