


# Precision

LABORATORIES





## RANGER DIPSLIDE INCUBATOR

**SKU:** INC-R45

Portable incubation of up to 45 dipslides. The Ranger Dipslide Incubator is designed for life in the field or in vehicles with its robust case resistant to oil and chemical environments. It can even be fit to an auxiliary battery so the incubator...

[Read More](#)



## CHALLENGER DIPSLIDE INCUBATOR

**SKU:** INC-C25

Portable incubation of up to 25 dipslides. The Challenger Dipslide Incubator case is easy to wipe clean and robust enough for most hostile work environments. The ergonomic location of the controls under the lid eliminates the risk of accidental...

[Read More](#)



## DISCOVERY DIPSLIDE INCUBATOR

**SKU:** INC-D8

Portable incubation of up to 8 dipslides. The Discovery Dipslide incubator is designed for life in the field or vehicles. The Temperature Control Technology is a fully solid state with no moving parts, and element switching is transistorized and...

[Read More](#)



## GNAT DIPSLIDE INCUBATOR

**SKU:** INC-G2

Portable incubation of 2 dipslides. The Gnat Dipslide Incubator is the easiest to operate. Just plug it in and within 30 minutes, it will heat up to the pre-set temperature and incubate until you turn it off. There are no additional buttons or...

[Read More](#)

